

Technical Manual

REBOUND® PCL

		3
EN		TECNICAL MANUAL 5
DE		GEBRAUCHSANWEISUNG 9
FR		NOTICE D'UTILISATION 14
ES		INSTRUCCIONES PARA EL USO 19
IT		ISTRUZIONI PER L'USO 23
NO		BRUKSANVISNING 28
DA		BRUGSANVISNING 32
SV		BRUKSANVISNING 36
EL		ΟΔΗΓΙΕΣ ΧΡΗΣΗΣ 40
FI		KÄYTTÖOHJEET 45
NL		GEbruIKSAANWIJZING 49
PT		INSTRUÇÕES DE UTILIZAÇÃO 53

ENGLISH

LIST OF PARTS (Figure 1)

- A. Proximal Thigh Strap
- B. Femoral Shell
- C. Distal Anterior Thigh Strap
- D. Hinge
- E. Condyle Pad
- F. Cable
- G. Cable Cover
- H. Shroud
- I. AMS Wrap
- J. Tibial Shell
- K. Distal Anterior Tibial Strap
- L. Upper Subshell
- M. Distal Posterior Thigh Strap
- N. Dynamic Tension System (DTS)
- O. Lower Subshell

Figure 7

- P. Spool
- Q. Shear Knob

Figure 9

- R. Pawl

INDICATION FOR USE

Rehabilitation of posterior cruciate ligament (PCL) injury, whether post-surgical reconstruction or during conservative, non-surgical treatment.

WARNINGS

- The brace must be prescribed and fit by a qualified healthcare provider.
- This brace is intended for post-surgical or post-injury rehabilitation. It is not intended to be used in moderate or high activity and/or impact sports. A qualified healthcare provider must approve the types of activities that the patient may engage in while using the brace.
- This knee brace is a supportive brace only and is not intended or guaranteed to prevent injury; treatment results may vary.
- The patient must not adjust the cable tension. Any adjustments to the cable tension must be made by the healthcare provider. The brace has been designed to be adjusted to the proper tension using the provided shear knobs. Adjustment using any other means can result in excessive tension and injury to the knee.
- In the event of accidental adjustment to the cable tension, the patient should immediately discontinue use of the brace and contact the healthcare provider.
- To avoid discomfort and potential impairment of circulation, the patient should avoid sitting on the frame for an extended period of time.
- If the patient experiences increased pain, swelling, sensation changes, excessive pressure, skin irritation or any adverse reactions while using this brace, immediately discontinue use of the brace and contact the healthcare provider.
- The brace must have an extension stop installed in the hinge during use.

- The patient should be instructed to keep the leg extended when sitting to increase comfort.

ENSURE THE END-USER AGREES TO THE FOLLOWING

1. Össur accepts no responsibility for any injury which might occur while using Össur knee bracing products.
2. With use of this type of device, it is possible to receive various types of skin irritation, including rashes and pressure points. If the brace causes discomfort or any unusual reactions, discontinue use and contact your health care provider immediately.
3. If you have any questions as to how this knee brace applies to your own specific knee condition, consult your health care provider.
4. This brace is not intended to be used in moderate or high activity and/or impact sports. If you do not agree to the above conditions, contact your brace provider before using this product.

INITIAL BRACE APPLICATION

1. Unfasten all straps and the AMS wrap.
Note: Brace comes installed with a 0° extension stop and a 105° flexion stop. Refer to the Range of Motion (ROM) adjustment section if change in ROM is required.
2. The leg must be in full extension when brace is applied. Place the hinge slightly above (approximately ¾”) the joint line so the brace can properly seat itself (**Figure 2**).
3. Beginning with the AMS wrap just below the knee, fasten the wrap snugly, as this will aid in brace suspension (**Figure 3**).
4. Stabilize the brace by loosely attaching the distal anterior thigh strap (C).
5. Next, snugly fasten the distal anterior tibial strap by inserting it through the slot in the lower subshell while centering the tibial shell on the Tibia (**Figure 4**).
6. Fasten the proximal thigh strap by inserting it through the slot in the upper subshell (**Figure 5**).
7. Fasten the anterior distal thigh strap and center the femoral shell so that the hinge is centered at midline (A/P) on the knee (**Figure 6**).
8. Fasten the distal posterior thigh strap to affix the hinge in its centered position at midline A/P.

SETTING INITIAL TENSION

9. Three different shear knobs are included:
 - a. White / three raised dots (standard load)
 - b. Grey / two raised dots (intermediate load)
 - c. Black / one raised dot (low load)**Note:** The standard load (white/three raised dots) is based on a review by Jansson et al KSSTA 2012
10. With patient in standing position and leg in full extension, place the desired shear knob in the diamond shaped hole on the spool. Slowly turn the shear knob in the clockwise direction until the shear knob breaks (**Figure 7**).
11. Discard the two pieces of the shear knob.
12. Place the DTS sticker on the spool to cover the diamond shaped hole. This is to prevent the patient and/or anyone other than the healthcare provider from adjusting the spool (**Figure 8**).
13. Provide instructions for use to the patient. Instruct the patient never to remove the DTS sticker and never to adjust the spool.

Note: All straps can be trimmed or adjusted to length. Remember, straps should not restrict circulation or cause cramping during use.

RE-TENSIONING DYNAMIC TENSIONING SYSTEM

1. Remove brace from leg if patient is still wearing the brace.
2. Unfasten the DTS system from the AMS wrap.
3. Flip the DTS system so the underside is exposed.
4. Insert a pen into the hole in the pawl and apply outward force in the direction of the arrow to disengage the pawl from the spool. While holding the pawl in the disengaged position, pull the cable out of the DTS, (this will unwind the cable and cause the spool to spin in the counter clockwise direction) until the DTS slides freely on the cable (**Figure 9**).
5. Center the DTS by pulling the DTS away from the brace and centering it between the uprights, (ensure that the cable is equal length on each side of the DTS) (**Figure 10**).
6. Place the desired shear knob into the diamond shaped hole on the spool and slowly turn the shear knob in the clockwise direction approximately a half turn to lock the DTS onto cable.
7. Flip the DTS back over and re-affix in its original position on the AMS wrap.

Note: The cables should be loose but of equal length on each side of the DTS.

8. Repeat steps 1–12 in the “*Initial Brace Application*” instructions to re-tension the DTS.

RANGE OF MOTION ADJUSTMENT (ROM)

Adjustment should be performed by a health care professional. The brace comes with 0° extension stops and 105° flexion stops attached to the hinges. If needed, ROM control can be adjusted using the extension and flexion stops included with the brace. Extension stops are labeled 0°, 5°, 10°, 20°, 30°, and 40°. Flexion stops are 0°, 15°, 25°, 35°, 45°, 55°, 65°, 75°, 90° and 105°.

ADJUSTING THE ROM

1. Remove the condyle pads from the medial and lateral hinges. Remove retaining screws. Using a screwdriver or pen, slide the 0° stop out of the hinge on both sides.
2. Re-apply the desired extension stop and gently secure it with the retaining screw (**DO NOT OVER TIGHTEN**).
3. Follow the same procedure to change the flexion stop.

CONDYLE PAD ADJUSTMENT

Attached to the hinges are the standard size pads. Pairs of thicker and thinner pads are included. The pads should lightly touch the medial and lateral condyles of the knee.

BRACE MAINTENANCE

The brace is designed for minimum maintenance and care. Follow these easy maintenance steps to keep the brace in good condition.

CLEANING

- Liners, straps and pads may be removed for cleaning.
- Do NOT machine wash.
- Do NOT place in dryer.
- Hand-wash using mild non-detergent soap.

- Rinse thoroughly.
- Air-dry.

If migration occurs, clean silicone liners with an alcohol wipe or damp cloth. If brace is used in salt water or chlorinated water, rinse well and air dry.

HINGE MAINTENANCE

- Do NOT disassemble the hinge
- Check regularly for foreign materials in hinge (i.e. dirt or grass)
- Carefully clean dirt from hinge using water or compressed air
- Lubricate with lithium grease

ACCESSORIES

The following accessories are available for your PCL brace at an additional charge. To order call your Össur Customer Service Department.

- PCL Flexion Stop kit
- PCL Extension Stop kit
- Replacement Liners
- Strap kit
- Condyle kit

LIMITED WARRANTY FOR END-USER

This warranty covers all defects in materials and workmanship of your new Rebound® PCL Knee Brace. This Warranty covers the frame and hinges for one year from date of sale. The warranty covers all other materials, including straps, sleeves, liners and pads for six months from the date of original purchase.

Össur will repair or replace (at its option) any defective part covered under warranty at no charge. Repairs necessitated by accident, negligence, misuse or abuse by the purchaser are not covered. Also, consequential and incidental damages are not recoverable under this warranty. These limitations may not apply as some states do not allow the exclusion or limitation of incidental or consequential damages. Only the original purchaser is entitled to claim service under this warranty.

To obtain service under this warranty call our Customer Service Department.

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- A. Proximaler Oberschenkel-Gurt
- B. Oberschenkel-Hülle
- C. Distal anteriorer Oberschenkel-Gurt
- D. Scharnier
- E. Femurkondylus-Auflage
- F. Kabel
- G. Kabelabdeckung
- H. Verkleidung
- I. AMS-Manschette
- J. Schienbein-Hülle
- K. Distal-anteriorer Schienbein-Gurt
- L. Obere Unterhülle
- M. Distal-posteriorer Oberschenkel-Gurt
- N. Dynamic Tension System (DTS) mit Spule
- O. Untere Unterhülle

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