

# Ovation Hybrid Night Splint

## APPLICATION INSTRUCTIONS

1. To apply, first unfasten all straps, **including tension strap**. Slide your foot into the splint from the rear (like a slipper) until it fits comfortable and snugly.




2. Apply the calf-strap by inserting ends through d-rings and tighten securely.



3. Apply foot-straps by pulling both straps up and over the foot to secure the splint snugly.



4. Flex the foot/ankle to desired position. Pull the tension strap downward and fasten.



Use position indicator to measure progress.



Approximately -10°



Approximately 90°



Approximately +10°

### Washing Instructions:

Wipe with damp cloth. Do not bleach. Line or towel dry only.

### CAUTION:

In the event of discomfort or numbness of any kind when using this product, consult your health care provider immediately.

The Hybrid Night Splint is not intended for continuous walking. Please consult your health care provider for specific instructions.



# Ovation Hybrid Night Splint

## APPLICATION INSTRUCTIONS

5. To apply, first unfasten all straps, **including tension strap**. Slide your foot into the splint from the rear (like a slipper) until it fits comfortable and snugly.




6. Apply the calf-strap by inserting ends through d-rings and tighten securely.



7. Apply foot-straps by pulling both straps up and over the foot to secure the splint snugly.



8. Flex the foot/ankle to desired position. Pull the tension strap downward and fasten.



Use position indicator to measure progress.



Approximately -10°



Approximately 90°



Approximately +10°

### Washing Instructions:

Wipe with damp cloth. Do not bleach. Line or towel dry only.

### CAUTION:

In the event of discomfort or numbness of any kind when using this product, consult your health care provider immediately.

The Hybrid Night Splint is not intended for continuous walking. Please consult your health care provider for specific instructions.



## Ovation Hybrid Night Splint - Accessory Strap (optional)

### APPLICATION SUGGESTION 1:



Figure 1

### APPLICATION SUGGESTION 2:




Figure 2

### APPLICATION SUGGESTION 3:




Figure 3

#### To use as additional support strap:

- Apply Ovation Hybrid Night Splint as instructed, then secure accessory strap to lower leg.

#### To provide extra support at toe strap:

- Apply Ovation Hybrid Night Splint as instructed, then secure accessory strap to toe strap.

#### To lengthen calf strap:

- Apply Ovation Hybrid Night Splint as instructed. Use accessory strap as extension for calf strap and cut excess strap not used.

Uses for accessory strap are not limited to these options. Patient can customize use of strap as needed or as directed by your health care provider.

### CAUTION:

In the event of discomfort or numbness of any kind when using this product, consult your health care provider immediately.



## Ovation Hybrid Night Splint - Accessory Strap (optional)

### APPLICATION SUGGESTION 1:



Figure 1

### APPLICATION SUGGESTION 2:




Figure 2

### APPLICATION SUGGESTION 3:



Figure 3

#### To use as additional support strap:

- Apply Ovation Hybrid Night Splint as instructed, then secure accessory strap to lower leg.

#### To provide extra support at toe strap:

- Apply Ovation Hybrid Night Splint as instructed, then secure accessory strap to toe strap.

#### To lengthen calf strap:

- Apply Ovation Hybrid Night Splint as instructed. Use accessory strap as extension for calf strap and cut excess strap not used.

Uses for accessory strap are not limited to these options. Patient can customize use of strap as needed or as directed by your health care provider.

### CAUTION:

In the event of discomfort or numbness of any kind when using this product, consult your health care provider immediately.

