

The Hybrid™ Night Splint (patent pending) combines the best elements of both PNS and DNS in treating Plantar Fasciitis with a unique, patient friendly design as adaptable as the patient's changing needs... Exclusive adjustable tensioning/positioning features result in increased patient compliance.

The comfort and convenience of a Dorsal Night Splint combined with the function of a Posterior Night Splint.

COMPONENTS



HYBRID™ MULTIPOSITIONAL NIGHT SPLINT



Hinge for ROM Foot component pivots up to 50 degrees for compound motion and adjustable positioning.



Position indicator Immediately shows the patient his/her position to help visualize goals and confirm progress. Printed in high contrast ink for low light situations, this component is designed with the patient in mind



Adjustable **tensioning strap (patent pending)** allows for individual custom adjustment. Each patient has a unique pain tolerance he/she experiences on a daily basis. This unique strap is the answer to quick, positional comfort and continued patient compliance.

MORE KEY FEATURES...



Wider foot component (4.25" width vs. typical 3") a full **42% wider** forefoot and toe area reduces toe cramping while comfortably increasing extended wear-time (resting)



Pivoting d-rings
Automatically adjusts to various calf shapes/profiles. Increases strap comfort and function.

Non-slip features A sticky, sure-grip feature on the soft weight bearing surface allows for safe ambulation and prevents slipping. **Foot grip** on the soft interfacing surface to prevent slipping.


APPLICATION

1. Unfasten ALL straps: 1 (calf strap), 2 (foot strap), and 3 (tension strap). Slide foot in like a slipper.

2. Apply calf strap (1) by inserting ends through d-rings and tighten securely.

3. Apply foot straps (2) by pulling both straps up and secure tightly over foot.

4. Flex foot/ankle to desired position using either option. Pull tension strap downward (to desired tension) and fasten.



SIZING

- 30014 S-M (Women 6 - 10; Men 5 - 9)
- 30016 L-XL (Women 10.5 - 15; Men 9.5 - 14)
- 30000S ALL

INDICATIONS FOR USE

- | | | |
|---------------------|------------|-------------------|
| Plantar fasciitis | Heel pain | Metatarsalgia |
| Achilles tendonitis | Heel spurs | Ankle contracture |
| Achilles tendinosis | Drop foot | Muscular strains |