



REBOUND® DUAL

Single source, complete collection



NEW!

Rebound® DUAL ST
(with step-through design)

REBOUND® DUAL

The Rebound DUAL is the versatile, comfortable, and functional support solution for knee ligament instability with, or without, osteoarthritis. Leveraging a light-weight profile and patented hi-tech ActiveGrip™ liners, the Rebound DUAL is suitable for low-to-high impact activities as well as sliding sports.



FEATURES AND BENEFITS

- 1 Sleek, strong frame**
Offers the strength of aircraft aluminum with a low-profile, cold-malleable frame that's easy to contour to patient's unique anatomy
- 2 Slide-to-size adjustability†**
Modify brace length from 12.5"-14" (32-35.5cm) or adjust varus / valgus angle
- 3 No-slip Activegrip™ calf liner**
Combines anti-gravity with breathability to reduce migration
- 4 Simple strapping system**
Easy-on, easy-off for the patient
- 5 Field-serviceable D-rings**
Technology for adjustments on the go
- 6 Comfortable control**
Doeskin thigh liner and cut-to-fit strap pads optimize patient comfort
- 7 Dynamic SmartDosing® technology‡**
Optional tool-free accessory for fine-tuned adjustment and offloading

† Slide-to-size Acculign frame available on Rebound DUAL only. Rebound DUAL ST is 14" (35.5 cm).

‡ Optional Smartdosing kit available for Rebound DUAL only, for ligament instabilities with osteoarthritis

REBOUND® DUAL ST

Evolving from the functional support and versatility of the Rebound DUAL in a new, step-through design, the Rebound DUAL ST was designed for improved ligament stability, exceptional comfort and fitting with ease.



CUSTOM COLORS

SHELL COLORS:

- Black (default)
- Electric Blue
- Metallic Silver
- White
- Red

LINER COLORS:

- Metallic Silver (default)
- Electric Blue
- Yellow
- Purple
- Black

Rebound DUAL Custom
with custom colors

ORDERING INFORMATION

Description	Size	Caliper ¹	Circumference ²	Part#
Rebound DUAL	XSmall	3"-3.5" (7.6-8.9 cm)	15"-16.5" (38.1-41.9 cm)	B-242XY0001
	Small	3.5"-4" (8.9-10.2 cm)	16.5"-18.5" (41.9-47 cm)	B-242XY0002
	Medium	4"-4.5" (10.2-11.4 cm)	18.5"-20.5" (47-52.1 cm)	B-242XY0003
or	Large	4.5"-5" (11.4-12.7 cm)	20.5"-22.5" (52.1-57.2 cm)	B-242XY0004
Rebound DUAL ST	X-Large	5"-5.5" (12.7-14 cm)	22.5"-24.5" (57.2-62.2 cm)	B-242XY0005
	XX-Large	5.5"-6" (14-15.2 cm)	24.5"-27.5" (62.2-69.9 cm)	B-242XY0006
	XXX-Large	6"-6.5" (15.2-16.5 cm)	27.5"-30" (69.9-76 cm)	B-242XY0007

Accessories	Rebound DUAL	Rebound DUAL ST	Measuring ³	Custom Part#
PCL Strap Kit XS-LG	B-701090057	B-701090065	SM, DMS, CCS, CAD, CAST	B-155XY0001
PCL Strap Kit XL-XXXL	B-701090058	B-701090066		
Smartdosing Kit RM/LL	B-701090060	N/A		
Smartdosing Kit RL/LM	B-701090061	N/A		

X = 5 (right) or 6 (left)

Y = 1 (Rebound DUAL) or 2 (Rebound DUAL ST)

¹ Caliper: M/L at joint width in a standard, weight-bearing position

² Thigh Circumference: 6" above mid-patella

³ SM = SmartMeasure, DMS = Digital Measuring System, CCS = Coordinate Cast System, CAD = Computer Aided Design

FOLLOW ÖSSUR ON



USA (800) 233-6263
WWW.OSSUR.COM

CANADA (800) 663-5982
WWW.OSSUR.CA

