



*Raising the Level of Patient Care*



# Adjustability Made Easy

The Vista Collar from Aspen improves patient care while saving time, money and storage space. With its innovative height adjustment technology, the Vista Collar is really six collars in one. The right size is always at hand, reducing storage and inventory costs, while improving patient care.



simplicity

## Eliminate Sizing

The Vista Collar eliminates the need for sizing methods of any kind. The innovative design of the Vista Collar allows the practitioner to simply apply the collar and then adjust to the proper height.

## Dial Height Adjustment

Getting the correct size for the patient is easy using the Vista Dial Height Adjustment™ System. Simply pull and turn the Dial to select one of six Vista height settings.

## Automatic Lock

Once the correct height has been determined, releasing the Dial automatically locks in the size.

patient care

## Effective Motion Restriction

Made by Aspen Medical Products, the Vista Collar provides safe and effective motion restriction<sup>1</sup> for all patients in need of extended wear collars.

## Improved Skin Care

The patented design of the Vista Collar disperses pressure, reducing the possibility of skin breakdown. Proven Aspen Pads<sup>2,3</sup> complete the system by wicking moisture away from the skin.

## Unparalleled Patient Access

The enormous tracheal aperture provides functional access for airway management and other essential procedures.

savings

## Inventory Management

By incorporating six sizes in one collar, the Vista Collar eliminates the cost of storing multiple sizes, frees up valuable inventory space, and eliminates SKUs.

## Reduce Waste

The Vista Collar reduces the waste associated with sizing errors and eliminates the need to change collars as patients become ambulatory.

## Save Time

Simple application and adjustment saves time in training and in use. "The Vista Collar improves your efficiency while improving your patient care."



Code L0174 Approved

Code L0172 Approved

1. Comparative Motion Allowed, Vista® vs. Aspen, Aspen Medical Products in-house testing, December 2006  
 2. Powers, J. A Multidisciplinary Approach to Occipital Pressure Ulcers Related to Cervical Collars. Journal of Nursing Quality, 1997, Vol. 2, No. 1  
 3. Absorption Comparison Test, Orange County Materials Test Laboratories, Anaheim CA, May 28,2002





# Adjustability for Support and Comfort

**T**he Vista Collar from Aspen improves patient care while saving time, money and storage space. With its innovative height adjustment technology, the Vista Collar is really six collars in one. The right size is always at hand, reducing storage and inventory costs, while improving patient care.




**T**he Vista TS from Aspen has multiple adjustment points to ensure a comfortable, effective fit for the vast variety of patient anatomies. Like the Vista Collar, the Vista TS is fully adjustable and only one size is required, reducing inventory costs and keeping patient care the main priority.



**Code L0172 Approved**

# Introducing the Vista® TS



simple

## One Adjustable Size

The Vista TS has several adjustment features that allow practitioners the ability to fit a wide variety of patient anatomies. Easy adjustments can be made on the Vista TS for a correct and comfortable fit for every patient, with only one stocking unit required.

## Multiple Adjustment Points

The angle adjustment mechanism provides a customized angle for each patient on the chest, while the height adjustment dial provides a custom length to each patient's torso.

## Easy Application

Metal free, the Vista TS is MRI compatible and requires no tools to fit or apply. Increased ease of application leads to increased patient compliance.



effective

## Adjustments for Anatomical Changes

Simply adjust the torso height dial to accommodate anatomical changes. This ensures that patients will be comfortable and supported while laying, sitting or standing.

## Stability and Skin Care

The Vista TS utilizes the Vista Collar so the patient benefits from the comfort and motion restriction of the collar, as well as the thoracic stabilizer.

## PRODUCT INFORMATION

| Vista Cervical Collar                                  |        | Vista TX Collar  |        | Vista TS Collar          |        |
|--|--------|------------------|--------|--------------------------|--------|
| Vista Cervical Collar                                  | 984000 | Vista TX Collar  | 984100 | Vista TS                 | 984500 |
| Vista Front Panel                                      | 984004 |                  |        | Vista TS Upgrade         | 984504 |
| Collar Sets <small>with extra replacement pads</small> |        | Replacement Pads |        | Back Panels              |        |
| Vista Cervical Collar                                  | 984002 | Vista Collar     | 984020 | Standard Back Panel      | 984006 |
| Vista TX Collar  | 984102 | Vista TX Collar  | 984120 | Standard Replacement Pad | 984026 |
|  |        | Vista TS         | 984520 | Large Back Panel         | 984008 |
|  |        |                  |        | Large Replacement Pad    | 984028 |



**Aspen Medical Products**  
 6481 Oak Canyon, Irvine, CA 92618-5202, USA  
 800-295-2776 • 949-681-0200  
[aspenmp.com](http://aspenmp.com)